



4021



1148



639



365



42



466



385



200



210



18



15



192



54

RUBBER GRIP  
PUSH HANDLES

CONTROLLER CAN BE FIXED  
ON BOTH SIDES

COMFORTABLE  
CUSHION

CLUTCH

SHOCKERS

ADJUSTABLE &  
DETACHABLE  
FOOTREST



REMOVABLE BATTERY BOX  
FOR EASY RECHARGING



STEEL WITH LIQUID  
PAINTING MATERIAL



MAG WHEELS WITH  
PNEUMATIC TYRES

### PARAMETERS

15-20KM  
DRIVING RANGE

DC 250W\*  
2 PCS

24 V 12AH  
(DOUBLE BATTERY)

PARKING  
BRAKE

ALUMINIUM  
ALLOY

MAX SPEED  
6KM/H

ELECTROMAGNETIC  
& MANUAL BRAKE

ATTENDANT  
CONTROLLER



# EVOX WC-102ME