

TM
AV Consumables



Say No to Specimen Recollections
Say Yes to Quality Backed by 30+ Years of
Diagnostic Lab Expertise


About AV Consumables

- AV Consumables was born out of the need to **avoid specimen recollections** for its parent organization, **Scientific Diagnostic Centre (SDC)**, Ahmedabad, India.
- Established in 1988, SDC is amongst India's first **NABL** certified pathology labs.
- **AV Consumable's leadership runs SDC's** pathology labs network, and has expertise in **Plastic engineering**
- AV's mission of "No Specimen Recollections for Labs" is ensured by **product innovation and commitment to quality**



Non-Vacuum Blood Collection Tubes

- Accurate Additives-Spraying by Robotic Technology
- Medical-grade butyl rubber caps
- Double-lock mechanism to ensure zero leakage



- Unbreakable and aesthetically pleasing
- Technician friendly – large labels make it easy to write
- Optimal packaging to reduce breakages in transportation

Additive

Volume

K3edta	2 & 3 ml
K2edta	2 & 3 ml
Clot Activator	4 ml
No Additive	4 ml
Sodium Fluoride + edta	2 ml
Lithium Heparin	3 ml
Gel + Clot activator	3.5 ml
Sodium Citrate	2 ml



Tubes Available in
Double or Rubber Cap



Thermocol Tray to Reduce
In-Transit Breakages

Ria Vials



- Made from medical-grade polypropylene (PP) and Polystyrene (PS)
- PS Ria tube is transparent, non-autoclavable and breakable
- PP Ria tube is Translucent, autoclavable and unbreakable
- Plastic push caps for longer storage

AV's Unique Packaging



- Neatly arranged; sealed in a pack of 100
- Minimizes transport cost and breakages

Serum Transport or Multi-purpose Vials



- False bottom allows tube to be loaded directly into all equipment
- Leak-proof and unbreakable
- Made from medical grade raw material



Flat and round bottom options available

Urine and Stool Containers

- Wide-mouth containers
- Made from medical-grade plastic
- Leak-proof screw cap mechanism
- Stool containers optimal for solid and semi-solid faeces
- Sterile and non-sterile options
- Container Sizes: Urine: 30, 40, 50 ml, Stool : 50 ml




Manufacturing Excellence

- AV has in-house manufacturing facilities in Ahmedabad, India
- A professional and well-trained team and advanced equipment ensures product excellence



Ferromatik Milacron machine produces high quality tubes



Nordson (USA) Spray Systems for precise spraying in tubes



CE and ISO 9001 Certified Operations



 **AV Consumables**

Phone: +91 9824 0911 05

Email: info@avconsumables.in

Website: www.avconsumables.in

Correspondence Address

3, Venunad,
Opp. White House, Panchwati,
Ahmedabad - 380006, Gujarat, India

Factory Address

69, Silicone Industrial Estate,
Opp. HOF Factory, Changodar,
Ahmedabad - 382213, Gujarat, India

