



MEDTECH[®]
Leader in Health & Home care!



OXYTEC-PRO[™]



OXYGEN CONCENTRATOR

Feature:

- Compact
- Light weight 16kg
- Low power- average 330 watts
- Quiet Operation
- Stainless Steel outlet
- Extensive Ventilation
- Can be used without air conditioning.
- 0-5 lpm upto 95% pure oxygen
- Valve and Motor American make
- Sieve from France

Useful:

- To Strengthen your heart, reducing the risk of heart attacks
- To Detoxify your Blood, Raise Energy Levels
- To Improve Strength & Build Endurance
- To Revitalize the skin
- To Prevent disease, relieve hangovers
- To Help fight off depression



MEDTECH[®]
Leader in Health & Home care!



1. Stainless Steel Oxygen outlet will not break during use.
2. Flow Indicator
3. Lockable flow Controller, in recessed cavity. Reduced chance of accidentally changing flow during use
4. ON/OFF Switch when ON, is inside the machine body, reduced chance of accidentally shutting off the machine
5. Reset switch on front panel
6. Front panel indicators for easy understanding of machine status.



Ventilation

The machine is very well ventilated. Air inlet for ventilation is in the front as well as the backside

The outlet for ventilation is also present in the front and the backside

Machine can be operated in non air conditioned room temperature environment also without fear of spoiling internal parts or sieve beds.

Oxygen is like gold

Oxygen is free but just like gold can be very hard to find. The most efficient and economical way to get enough oxygen is in the way we breathe.

It is advised to learn to breathe better, eat more live enzyme, nutrient and oxygen-rich foods and exercise moderately without excessive straining, gasping or breath heaving.

MEDTECH LIFE PVT. LTD.

(Formerly Nulife Global Medical Devices Pvt. Ltd.)

B6, Byculla Service Industries,
D.K. Marg, Byculla, Mumbai 400 027, India.
T: +91 22 23748372/73 / +9122 6657 8989
E: info@medtechlife.com
W : www.medtechlife.com



An ISO 9001: 2008
Certified Company