



EVOX

WHEEL CHAIR SERIES

Contributing To Healthy Life

EVOX WC-107

POWER WHEEL CHAIR



Features: Foldable directly
Removable battery box for recharging
Controller can be exchanged from right side to the left

Parameters: ● Driving range : 15-20km

- Motor: DC 250W*2Pcs
- Battery: 24V 20Ah (Lithium battery)
- Parking Brake: **Electromagnetic brake**
- Frame material: Aluminium alloy
- Frame Color: ○ silver
- Cushion: ● Black
- Max Speed: 6km/h



27/28	980	600	950	445	460	460	460	N/A	8	12	100