

# EVOX WC-104



- Features:** Extra function on
- .Electric elevating footrest: adjustable angle 60
  - .Electric reclining backrest: adjustable angle 45



## Parameters: ● Driving range : 15-20km

- Motor: DC 250W\*2Pcs
- Battery: 24V20Ah
- Parking Brake: **manual brake**
- Frame material: Steel with liquid painting
- Frame Color: ○ silver,
- Cushion: ● Black
- Max Speed: 6km/h



N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)	backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)