



EVOX

Contributing To Healthy Life

WHEEL CHAIR SERIES

EVOX WC-107

POWER WHEELCHAIR



Features:

- Easy Folding without removing any accessory
- Very Light in weight
- Removable Battery box for easy recharging
- Controller can be fixed on both sides
- Attendant Controller (Optional)

Parameters: ● Driving range : 15-20km

- Motor: DC 250W*2Pcs
- Battery: 24V 20Ah (Lithium battery)
- Parking Brake: **Electromagnetic brake**
- Frame material: Alumminium alloy
- Frame Color: ○ silver
- Cushion: ● Black
- Max Speed: 6km/h



N.W. (kg)	overall length (mm)	overall width (mm)	overall height (mm)	seat width (mm)	seat depth (mm)	backrest height (mm)	seat height (mm)	folded width (mm)	front wheel (inch)	rear wheel (inch)	loading capacity (kg)
27/28	980	600	950	445	460	460	460	N/A	8	12	100

6238810220



- Features:**
- .Foldable Dual Operation : Manual and Electric Mode
 - .Removable Battery box for easy recharging
 - .Flip-up armrest, adjustable and detachable footrest
 - .Controller can be fixed on both sides
 - .Handle Height Adjustment (Optional)
 - .Seat Width 18/20 (Optional)
 - .Attendant Controller (Optional)



Parameters: ● Driving range : 15-20km

- Motor: DC 250W*2Pcs
- Battery: 24V 12Ah (Double Battery)
- Parking Brake: manual brake
- Frame material: steel with liquid painting
- Frame Color: ○ silver
- Cushion: ● Black & red,
- Max Speed: 6km/h



N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)	backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
40/46	1090	685	930	450	435	395	500	420	8	22	100



- Features:**
- Foldable Dual Operation : Manual and Electric Mode
 - Removable Battery box for easy recharging
 - Flip-up armrest, adjustable and detachable footrest
 - Controller can be fixed on both sides
 - Attendant Controller (Optional)
 - With Electromagnetic and manual breaks



Parameters: ● Driving range : 20-25km

- Motor: DC 250W*2Pcs spur gear motor
- Battery: 24V 12Ah (Double battery)
- Parking Brake: Manual
- Frame material: Steel with liquid painting
- Frame Color: ○ Silver
- Cushion: ● Black
- Max Speed: 6km/h



N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)	backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
45/51	1140	630	980	450	430	395	500	370	10	16	100



Features: Light weight for easy carrying
 Foldable Dual Operation : Manual and Electric Mode
 Removable Battery box for easy recharging
 Flip-up armrest, adjustable and detachable footrest
 Controller can be fixed on both sides
 Attendant Controller (Optional)



Parameters: ● Driving range : 15-20km
 ● Motor: DC 250W*2Pcs
 ● Battery: 24V 12Ah (Lithium Battery)
 ● Parking Brake: manual brake
 ● Frame material: aluminium alloy
 ● Frame Color: ● Blue
 ● Cushion: ● Blue mess
 ● Max Speed: 6km/h

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)	backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
36/41	1030	685	955	450	435	395	500	415	8	22	100

Features:
 Reclining Electric WheelChair
 Reclining backrest and footrest position can be achieved with remote control
 High performance Taiwan Imported Motor with inbuilt Electromagnetic breaks
 Removable Battery box for easy recharging
 Flip-up armrest, adjustable and detachable footrest
 Controller can be fixed on both sides
 Attendant Controller (Optional)



Parameters: ● Driving range : 15-20km
 ● Motor: DC 250W*2Pcs (Taiwan)
 ● Battery: 24V20Ah(Double Battery)
 ● Parking Brake: Manual & Electromagnetic brake
 ● Frame material: Steel with liquid painting
 ● Frame Color: ○ silver,
 ● Cushion: ● Black
 ● Max Speed: 6km/h

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)	backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
57	1215	620	940	450	435	395	500	370	8	12	120



Features:

- Foldable
- Removable Battery box for easy recharging
- Flip-up armrest, adjustable and detachable footrest
- Controller can be fixed on both sides
- Electromagnetic Break (Optional)
- Attendant Controller (Optional)



Parameters: ● Driving range : 15-20km

● Motor: DC 250W*2Pcs

● Battery: 24V20Ah (Double Battery)

● Parking Brake: manual & electromagnetic brake

● Frame material: steel with liquid painting

● Frame Color: ○ silver

● Cushion: ● Black & red,

● Max Speed: 6km/h

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)	backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
41/46 42/47(E)	1140	640	935	450	435	395	500	370 470(E)	8	12	100

Attendant Controller



Features: Attendant controller is a dual control device designed for attendant/caregivers to help the user while walking behind the power wheelchair.

The attendant controller specially designed for exclusively for the PWC use only.

The attendant controller is attached close to hand grips

The attendant controller can be connected with joystick Knob of every EVOX-WC models

Provide the best possible care for attendant/caregivers while achieving normal driving speed.

Very easy to handle & operate